

# COMPARISON CHART

Brand Name	Liquiband®	SwiftSet®	SurgiSeal®	SurgiSeal® Twist	Skin Affix®	Liquiband® Octyl	McKesson Liquiband® Exceed	Dermabond Advanced®	Dermabond® Mini	exofin® HIGH VISCOSITY TISSUE ADHESIVE
<b>FDA 510k Name</b>	Liquiband Flow Control	Liquiband Flow Control	SurgiSeal	SurgiSeal Stylus	Barle Tissue Adhesive 2	Barle Tissue Adhesive 2	Barle Tissue Adhesive 2	Dermabond NX	Dermabond NX	Exofin HV Tissue Adhesive
<b>Premarket Submission #</b>	K122446	K122446	K082993	K130474	K132243	K132243	K132243	K100423	K100423	K152476
<b>Manufacturer</b>	Advanced Medical Solutions	Advanced Medical Solutions	Adhezion Biomedical	Adhezion Biomedical	Advanced Medical Solutions	Advanced Medical Solutions	Advanced Medical Solutions	Ethicon, Inc.	Ethicon, Inc.	Chemence Medical, Inc.
<b>Distributor</b>	Cardinal Healthcare	Covidien Ltd	Medline Industries	Medline Industries	Medline Industries	Cardinal Healthcare	McKesson	Ethicon, Inc.	Ethicon, Inc.	Chemence Medical, Inc.
<b>Distributed by Manufacturer</b>	—	—	—	—	—	—	—	✓	✓	✓
<b>Composition/Active Monomer</b>	n-butyl	n-butyl	2-Octyl	2-Octyl	2-Octyl	2-Octyl	2-Octyl	2-Octyl	2-Octyl	2-Octyl
<b>Size (Fill)</b>	.7ml	.8ml	.36ml	.5ml	.4ml	.8ml	.8ml	.7ml	.4ml	1.0ml
<b>Viscosity**</b>	~ 5 cPs	~ 5 cPs	~ 13 cPs	~ 13cPs	~ 59 cPs	~ 59 cPs	~ 59 cPs	~ 266 cPs	~ 266 cPs	~ 450 cPs
<b>Polymerization Time &lt; 1 Min**</b>	✓	✓	—	—	—	—	—	—	—	✓
<b>Flexible on Skin</b>	—	—	✓	✓	✓	✓	✓	✓	✓	✓
<b>Waterproof</b>	*	*	✓	✓	✓	✓	✓	✓	✓	✓
<b>Microbial Barrier</b>	*	*	✓	✓	✓	✓	✓	✓	✓	✓
<b>Transparent Applicator</b>	—	—	—	—	—	—	—	✓	✓	✓
<b>Flexible Applicator Tip</b>	—	—	—	—	—	—	—	—	—	✓

\* n-butyl monomer is prone to crack and break on skin thus compromising waterproof and microbial barrier

\*\* Data on file

